


 <p> <b>UINVERSITAS ISLAM NEGERI</b>  <b>SYARIF HIDAYATULLAH</b>  <b>FAKULTAS SAINS DAN TEKNOLOGI</b>          Jl. Ir. H. Juanda No. 95 Ciputat 15412 Indonesia       </p>	<b>STANDAR OPERASIONAL PROSEDUR (SOP)</b>	No. Dok : _____
		Tgl. Terbit : 12 Juli 2025
		No. Revisi : _____
		Hal : 1/2
<b>GLASSWARE CLEANING WORK INSTRUCTIONS</b>		

1. Ensure that you are wearing appropriate safety gear, such as gloves and eye protection.
2. Prepare a clean and disturbance-free work area.
3. Gather all glassware to be cleaned.
4. Separate glassware contaminated with hazardous chemicals.
5. Use distilled water to rinse the glassware.
6. If scrubbing is necessary, use a brush with soapy water.
7. Rinse with tap water, followed by a rinse with distilled water.
8. Use ethanol to clean contaminants soluble in ethanol from the glassware.
9. Follow with a rinse of distilled water.
10. Rinse with any other solvent needed, followed by ethanol and distilled water.
11. After cleaning, carefully dry the glassware.
12. Store the glassware in a clean and appropriate location.

 <p> <b>UIN</b>  <b>UNIVERSITAS ISLAM NEGERI</b>  <b>SYARIF HIDAYATULLAH</b>  <b>FAKULTAS SAINS DAN TEKNOLOGI</b>  <small>Jl. Ir. H. Juanda No. 95 Ciputat 15412 Indonesia</small> </p>	<b>STANDAR OPERASIONAL PROSEDUR (SOP)</b>	No. Dok : _____
		Tgl. Terbit : 12 Juli 2025
		No. Revisi : _____
		Hal : 2/2
<b>GLASSWARE CLEANING WORK INSTRUCTIONS</b>		

Disahkan oleh:	Disiapkan oleh:
Kabag TU	Mahasiswa
	
<b>Dra. Faojah, M.A.</b>	<b>Anjani Khoerul Fadilah</b>
NIP: 196803031996032001	NIM: 11220860000111